


RV INVESTIGATOR

HYDROCHEMISTRY DATA PROCESSING REPORT

Voyage:	IN2023_V03
Chief Scientist	Elizabeth Shadwick (CSIRO / UTAS)
Voyage title:	SOTS: Southern Ocean Time Series automated moorings for climate and carbon cycle studies southwest of Tasmania
Report compiled by:	Harris Anderson & Narendra Pati



Contents

1	Executive Summary	4
1.1	Objectives.....	4
1.2	General Hydrochemistry Information.....	4
2	Itinerary	6
3	Key personnel list	6
4	Summary.....	7
4.1	Sample Type and Number Assayed.....	7
4.1.1	CTD samples (Conductivity, Temperature, Density)	7
4.1.2	Thermosalinograph (TSG) samples	7
4.2	Data Processing Overview	8
4.2.1	Conventional hydrology data.....	8
5	Salinity	9
5.1	Salinity Measurement Parameters	9
5.2	Salinity Method.....	9
5.3	CTD Salinity vs Bottle Salinity Plot	10
6	Dissolved Oxygen.....	11
6.1	Dissolved Oxygen Measurement Parameters.....	11
6.2	Dissolved Oxygen Method	11
6.3	CTD Dissolved Oxygen vs Bottle Dissolved Oxygen Plot	12
6.4	Dissolved Oxygen Instrument titrant: thiosulphate normality and blank correction.....	13
7	Nutrients.....	15
7.1	Nutrient Measurement Parameters	15
7.2	Nutrient Methods	16
7.3	HyPro Processing Summary for Nutrients	16
7.4	Accuracy - Reference Material for Nutrient in Seawater (RMNS)	18
7.5	Nutrient plots of RMNS.....	18
7.6	Measurement Uncertainty.....	21
7.7	Method Detection Limit for Nutrients.....	21
7.8	Sampling Precision	22
7.9	Redfield Ratio Plot (14.0) for CTD Deployments.....	23
7.10	Temperature & Humidity Change over Nutrient Analyses	23
8	Appendix.....	25
8.1	Salinity: Reference material used	25

8.2	Nutrients: RMNS results for each CTD Deployment.....	25
8.2.1	Lot CP	25
8.3	Measured MDL for each CTD deployment.....	26
8.3.1	AA3HR Auto analyser	26
8.4	Missing or Suspect Salinity Data	26
8.5	Missing or Suspect Dissolved Oxygen Data	26
8.6	Missing or Suspect Nutrient Data	27
8.7	Data Quality Flag Key	27
8.8	GO-SHIP Specifications.....	27
8.8.1	Salinity.....	27
8.8.2	Dissolved Oxygen	28
8.8.3	Si(OH) ₄	28
8.8.4	PO ₄	28
8.8.5	NO ₃	28
8.8.6	Notes	28
9	References	29

1 Executive Summary

1.1 Objectives

The primary objective for in2023_v03 was to first deploy a new set of SOTS moorings (SOFS-12 and SAZ-25) and then recover the existing SOTS moorings (SOFS-11 and SAZ-24). Each of the SOTS moorings delivers to specific aspects of the atmosphere-ocean exchanges:

- the SAZ sediment trap mooring collects samples to quantify the transfer of carbon and other nutrients to the ocean interior by sinking particles and investigate their ecological controls.
- the Southern Ocean Flux Station (SOFS) mooring measures meteorological and ocean properties important to air-sea exchanges, ocean stratification, waves, currents and biological productivity and ecosystem structure. Water samples are collected for more detailed nutrient and plankton investigations after recovery.

Ancillary work obtained supporting information on atmospheric and oceanographic conditions using CTD casts, underway measurements, Continuous Plankton Recorder, and autonomous glider, and potentially casts of a bio-optical sensor package.

Five total CTDs deployed, including four at the SOTS/SOFS mooring site. CTDs sampled and analysed by Hydrochemistry team.

1.2 General Hydrochemistry Information

Water samples collected during the voyage were analysed in the ship's hydrochemistry laboratory for nutrients, dissolved oxygen, and salinity. Overall data collected was of high quality. No significant sample collection, analysis, or data processing issues were encountered.

Five nutrients were determined: silicate, phosphate, nitrate + nitrite, nitrite and ammonium using AA3HR autoanalyser. Certified reference materials for nutrients in seawater (RMNS) were within 3% of their certified values. Missing and suspect hydrology samples are listed in Appendix section.

Please cite the following manuscript when reporting or publishing data for silicate, phosphate, nitrate+nitrite (NO_x) and nitrite:

Rees, C., L. Pender, K. Sherrin, C. Schwanger, P. Hughes, S. Tibben, A. Marouchos, and M. Rayner. (2018) "Methods for reproducible shipboard SFA nutrient measurement using RMNS and automated data processing."

Limnol. Oceanogr: Methods, 17(1): pp. 25-41.

doi:10.1002/lom3.10294

If publishing ammonium data, please cite the following:

Rees, C., Janssens, J., Sherrin, K., Hughes, P., Tibben, S., McMahon, M., McDonald, J., Camac, A., Schwanger, C. and Marouchos, A., (2021) *“Method for Reproducible Shipboard Segmented Flow Analysis Ammonium Measurement Using an In-House Reference Material for Quality Control.”*

Frontiers in Marine Science, 8.

doi:10.3389/fmars.2021.581901

Final hydrology data, analytical methods, related log sheets and processing notes can be obtained from the CSIRO data centre.

For Data, contact: [NCMI DataLibrarians@csiro.au](mailto:NCMI_DataLibrarians@csiro.au)

2 Itinerary

Table 1: Voyage itinerary

	Depart	Arrive
Port	Hobart	Hobart
Date	12/05/2023	25/05/2023
Time	1100	0830




Figure 1. Voyage track

3 Key personnel list

Table 2: Key Personnel list

Name	Role	Organisation
Elizabeth Shadwick	Chief Scientist	CSIRO/UTAS
Margot Hind	Voyage Manager	CSIRO
Rod Palmer	Alternate Voyage Manager	CSIRO
Harris Anderson	Hydrochemist	CSIRO
Narendra Pati	Hydrochemist	CSIRO

4 Summary

4.1 Sample Type and Number Assayed

Table 3: Sample Type and Number Assayed

Analysis	Samples Assayed	Type
Salinity	77	CTD
	12	TSG
Dissolved Oxygen	77	CTD
Nutrients (AA3)	74	CTD
	23*	exp

*Note that some experimental samples were re-analysed for Phosphate after the voyage (05/05/2023)

4.1.1 CTD samples (Conductivity, Temperature, Density)

- Taken from the 12L Ocean Test Equipment bottles on the CTD rosette that is deployed at depth for water collection.
- A total of 5 CTD deployments were sampled by
 - Hydrochemistry: Harris Anderson and Narendra Pati
 - Science party: Patrick Duke

4.1.2 Thermosalinograph (TSG) samples

- Taken from the underway instrument clean seawater line supplying the pCO₂ instrument in the underway laboratory.
- TSG samples collected by hydrochemistry. Results emailed to Vito Dirita (CSIRO) at the completion of the voyage.

Refer to voyage EVERLog for UWY and TSG sample information.

4.2 Data Processing Overview

4.2.1 Conventional hydrology data

The sample meta-data, measured bottle salinity results, dissolved oxygen assay results and the nutrient assay raw data are processed by the CSIRO program HyPro. The final output is the hydrology data set. An overview of this process is illustrated below (fig.2).




Figure 2. Conventional Hydrology Data Processing Flow Diagram.

5 Salinity

5.1 Salinity Measurement Parameters

Table 4: Salinity Measurement Parameters

Details	
HyPro Version	5.7
Instruments	Guildline Autosal Laboratory Salinometer 8400(B) – SN 72088. Bath temperature 24.0°C
Software	Ocean Scientific International Ltd (OSIL) Data Logger ver 1.2
Hydrochemistry Methods	Sampling: WI_Sal_002 Analysis : SOP 006
Accuracy	± 0.001 practical salinity units
Reference Material	OSIL IAPSO – Batch P166, use by 06/04/2025, $K_{15} = 0.99987$
Sample Container	200 mL volume OSIL bottles made of type II glass (clear) with disposable plastic insert and plastic screw cap.
Sample Storage	Stored in salinometer lab for minimum of 8 hrs before measurement.
Lab Temperature	Mean 20.9°C SD 0.9°C
Analysts	Harris Anderson
Comments	See DAP report for CTD calibration details.

5.2 Salinity Method

Salinity samples were measured on a Guildline Autosal 8400B instrument operated in accordance with its technical manual. The measured value is recorded with an OSIL data logger.

Practical salinity (S) is defined in terms of the ratio (K_{15}) of the electrical conductivity measured at 15°C 1atm of seawater to that of a potassium chloride (KCl) solution of mass fraction 32.4356×10^{-3} .

Before each lot of sample measurements, the Autosal is calibrated with standard seawater (OSIL, IAPSO) of known K_{15} ratio. A new bottle of OSIL standard is used for each calibration. The frequency of calibration is at least one per set of samples per CTD deployment.

Method: The salinity sample is collected in a 200ml OSIL bottle. The bottle is rinsed then filled from the bottom, via a polytetrafluoroethylene (PTFE) straw, till overflowing. The bottle is removed from the straw and the sample is decanted to allow a headspace of approximately 25cm³. A dry plastic insert is fitted, the bottle inverted and rinsed with water then capped and stored cap-down until measured. To measure, the Autosal cell is flushed three times with the sample and then measured after the fourth and fifth flush. The OSIL data logger software captures the conductivity ratio and calculates the practical salinity.

The output from the data logger is imported into HyPro and collated with the CTD deployment meta-data.

5.3 CTD Salinity vs Bottle Salinity Plot

For this voyage, the difference between the unprocessed (uncorrected) CTD value and the measured bottle value is generally less than 0.002 PSU. The larger differences are for shallow samples across the sudden changes in the thermohaline profile.

The unprocessed CTD values are adjusted (corrected) by DAP using the bottle results. The corrected values are not reported in the hydrology set. Please contact the [NCMI DataLibrarians@csiro.au](mailto:NCMI_DataLibrarians@csiro.au) for corrected CTD data.




Figure 4. CTD Salinity - Bottle Salinity vs CTD deployment plot. The data quality is coded by colour and delineated by a dot for the bottle salinity and a circle for the CTD salinity. Green = GOOD. Black = UNPROCESSED. Units: PSU (dimensionless).

6 Dissolved Oxygen

6.1 Dissolved Oxygen Measurement Parameters

Table 5: Dissolved oxygen measurement parameters.

Details	
HyPro Version	5.7
Instrument	Scripps Automated Photometric Oxygen System (SIO)
Software	LVO2 ver 2.36 Scripps Institution of Oceanography (SIO)
Hydrochemistry Methods	Sampling: WI_DO_001 Analysis: SOP 005
Accuracy	$\pm 0.5 \mu\text{mol L}^{-1}$
Analysts	Narendra Pati
Lab Temperature ($\pm 1^\circ\text{C}$)	Mean 17.7°C SD 0.8°C
Sample Container type	140 mL glass iodine determination flasks with glass stopper.
Sample Storage	Samples stored in the hydrochemistry lab until analysis. All samples were analysed within ~48 hrs
Comments	See DAP report for CTD calibration details.

6.2 Dissolved Oxygen Method

SIO method used. The method is based on the whole bottle modified Winkler titration of Carpenter (1965) plus modifications by Culberson *et al* (1991).

Method: The sample is collected in an iodine determination flask of known volume. 1mL of manganese (II) chloride solution followed by 1 mL of alkaline iodide solution is added to the sample, the flask stoppered and inverted a minimum of 15 times. The dissolved oxygen oxidizes an equivalent amount of Mn (II) to Mn (IV) which precipitates. Just before titration, the sample is acidified, Mn (IV) is reduced to the divalent state liberating iodine. The iodine is titrated with a standardised thiosulphate solution using a Metrohm 665 Dosimat fitted with a 1 mL burette. The endpoint is determined by measuring the decrease in the UV absorption 365 nm.

The thiosulphate solution is standardised by with a 10 mL aliquot of potassium iodate primary standard. A blank correction is also determined from the difference between two titres of consecutive additions of 1 mL aliquots of potassium iodate to the same blank sample. The standardisation is done at least once per 12-hour shift, when samples are being assayed.

The output from the SIO instrument software is imported into HyPro and collated with the CTD deployment meta-data.

6.3 CTD Dissolved Oxygen vs Bottle Dissolved Oxygen Plot

For this voyage, the difference between the unprocessed CTD value and the measured bottle value is generally less than 15 $\mu\text{mol L}^{-1}$. The larger differences are for shallow samples across the sudden changes in the dissolved oxygen profile.

The unprocessed CTD values are adjusted (corrected) by DAP using the bottle results. The corrected values are not reported in the hydrology set. Please contact the [NCMI DataLibrarians@csiro.au](mailto:NCMI_DataLibrarians@csiro.au) for corrected CTD data.




Figure 5. CTD Dissolved Oxygen - Bottle Dissolved Oxygen vs Deployment Plot. The data quality is coded by colour and delineated by a dot for the bottle DO and a circle for the CTD DO. Green = GOOD. Blue = SUSPECT. Red = BAD. Black = UNPROCESSED. Units: $\mu\text{mol L}^{-1}$. *Note: Bad oxygen bottle data is listed in appendix 8.5.

6.4 Dissolved Oxygen Instrument titrant: thiosulphate normality and blank correction.

The variance in thiosulphate concentration is within our QC parameter of less than 0.0005 N between standardisations. The blank correction is used in the calculation of the thiosulphate normality and is due to oxidisable species in the MQ water that is added to the KIO_3 aliquot before the titration.

For thiosulphate normality plots, the red lines indicate ± 0.0005 N either side of the mean titrant (thiosulfate) concentration. For blank plot, red lines indicate acceptable variation either side of the mean blank concentration. The titrant should not vary more than 0.0005 N between analyses.




Figure 6. Thiosulphate standardisation and blank correction plots.

7 Nutrients

7.1 Nutrient Measurement Parameters

Table 6: Nutrient measurement parameters analysed with Seal AA3HR segmented flow analyser. All instrument parameters, reagent batches and instrument events are logged for each analysis run. This information is available on request.

Details					
Instrument	Seal AA3HR segmented flow analyser				
HyPro version	5.7				
Operating Software	AACE 7.10				
Hydrochemistry Sampling Method	WI_Nut_001				
Hydrochemistry analysis method	SOP001	SOP002	SOP003	SOP003	SOP004
Nutrients	Silicate	Phosphate	Nitrate + Nitrite	Nitrite	Ammonia
Top concentration ($\mu\text{mol L}^{-1}$)	140.0	3.0	42.0	1.4	2.0
Method detection limit ($\mu\text{mol L}^{-1}$)	0.2	0.02	0.02	0.02	0.02
Reference Material	KANSO RMNS lot CP				
Sample Container	CTD: 50 mL HDPE with screw cap lids. Reused after acid wash with 10% HCl solution. UWY: 12 mL PP tubes with screw cap lids.				
Sample Storage	< 4 hours at room temperature after collection or < 12 hours at 4°C after collection				
Sample preparation	Assayed as neat. No filtration.				
Lab Temperature (°C)	Mean 17.7°C SD 0.8°C				
Analysts	Harris Anderson and Narendra Pati				
Comments	N/A				

7.2 Nutrient Methods

Nutrient samples are assayed on a Seal AA3HR segmented flow auto-analyser fitted with 1cm flow-cells for colorimetric measurements and a JASCO FP2020 fluorescence instrument as the ammonium detector.

Silicate (SOP001): colourimetric, molybdenum blue method. Based on Armstrong et al. (1967). Silicate in seawater is reacted with acidified ammonium molybdate to produce silicomolybdic acid. Tartaric acid is added to remove the phosphate molybdic acid interference. Tin (II) chloride is then added to reduce the silicomolybdic acid to silicomolybdous acid and its absorbance is measured at 660nm.

Phosphate (SOP002): colourimetric, molybdenum blue method. Based on Murphy and Riley (1962) with modifications from the NIOZ-SGNOS¹ Practical Workshop 2012 optimizing the antimony catalyst/phosphate ratio and the reduction of silicate interferences by pH. Phosphate in seawater forms a phosphomolybdenum complex with acidified ammonium molybdate. It is then reduced by ascorbic acid and its absorbance is measured at 880nm.

Nitrate (SOP003): colourimetric, Cu-Cd reduction – naphthylenediamine method. Based on Wood et.al (1967). Nitrate is reduced to nitrite by first adding an ammonium chloride buffer then sending it through a copper - cadmium column. Sulphanilamide is added under acidic conditions to form a diazo compound. This compound is coupled with 1-N-naphthly-ethylenediamine di-hydrochloride to produce a reddish purple azo complex and its absorbance is measured at 520 nm.

Nitrite (SOP003): colourimetric, naphthylenediamine method. As per nitrate method without the copper cadmium reduction column and buffer.

Ammonium (SOP004): fluorescence, ortho-phtaldialdehyde method. Based on K rouel and Aminot (1997). Ammonium reacted with ortho-phtaldialdehyde and sulphite at a pH of 9.0-9.5 to produce an intensely fluorescent product. Its emission is measured at 460nm after excitation at 370nm.

SOP methods can be obtained from the CSIRO Oceans and Atmosphere Hydrochemistry Group.

¹ Royal Netherlands Institute for Sea Research – Study Group on Nutrient Standards.

7.3 HyPro Processing Summary for Nutrients

After a run, the raw absorbance/ fluorescence data is exported from the instrument and processed by HyPro. For each analyte, HyPro re-creates the peak traces, defines the region on the peak's plateau (peak window) used to determine the peak heights, constructs the calibration curve, applies corrections for carry-over, baseline, and sensitive drifts then, derives the nutrient concentrations for each sample. The corrections are quantified using dedicated solutions included in every run.

HyPro uses criteria to identify suspect calibration points, noisy peaks, method detection limits that are above the nominal limit and duplicate sample results that do not match.

Suspect calibration points are weighted less when fitting the calibration curve. The cut-off limits for good calibration data are:

- $\pm 0.5\%$ of the concentration of the top standard for silicate and nitrate+nitrite (as per WOCE¹).
- $0.02 \mu\text{mol L}^{-1}$ for phosphate, nitrite, and ammonium.

HyPro classifies the quality of data as good, suspect, or bad and flags accordingly. The Flag key is in Appendix 8.7. Missing or suspect nutrient data is tabulated in Appendix 8.6.

¹ World Ocean Circulation Experiment

Table 8: HyPro 5.7 Processing Parameters. All instrument parameters and reagent batches and operation events are logged for each analysis run. This information is available on request.

Result Details	Silicate	Phosphate	Nitrate + Nitrite (NOx)	Nitrite	Ammonia
Data Reported as	$\mu\text{mol L}^{-1}$	$\mu\text{mol L}^{-1}$	$\mu\text{mol L}^{-1}$	$\mu\text{mol L}^{-1}$	$\mu\text{mol L}^{-1}$
Calibration Curve fit	Linear	Linear	Linear	Quadratic	Quadratic
# of points in Calibration	6	6	6	6	6
Forced through zero	N	N	N	N	N
Matrix correction	N	N	N	N	N
Blank correction	N	N	N	N	N
Peak window defined by	HyPro	HyPro	HyPro	HyPro	HyPro
Carryover correction (HyPro)	Y	Y	Y	Y	Y
Baseline drift correction (HyPro)	Y	Y	Y	Y	Y
Sensitivity drift correction (HyPro)	Y	Y	Y	Y	Y
Data Adj for RMNS variance.	N	N	N	N	N
Medium of Standards	Low nutrient seawater (LNSW, bulk on PW1 wharf, CSIRO Hobart) collected in June 2021. Sub-lot passed through a 10-micron filter (filtered on 14/03/2023) and stored in 20 L carboys in the clean dry laboratory at 22°C.				
Medium of Baseline	18.2 Ω water. Dispensed from the Milli Q IQ 7010 system.				
Duplicate samples	CTD: Niskin fired at the greatest depth were analysed in duplicate. Single samples were analysed for remaining depths.				
Comments	The reported data is not corrected to the RMNS. Per deployment RMNS data tabulated in appendix 8.2.				

Table 9: HyPro 5.7 Processing Parameters. All instrument parameters and reagent batches and operation events are logged for each analysis run. This information is available on request.

7.4 Accuracy - Reference Material for Nutrient in Seawater (RMNS)

Descriptive statistics are used to ascertain the accuracy and precision of the analysis from the repetitive measurement of the RMNS for silicate, phosphate, NO_x, and nitrite in seawater.

For in2023_v03, the certified reference material results lot CP for NO_x and nitrite are within 1% of their respective certified mean concentration, whereas silicate and phosphate are within 2% of their respective certified mean concentration.

The GO-SHIP criteria (Hyde *et al.*, 2010), appendix 8.8, specifies using 1-3 % of full scale (depending on the nutrient) as acceptable limits of accuracy.

KANSO publishes the RMNS nutrient values in $\mu\text{mol kg}^{-1}$. These are converted to $\mu\text{mol L}^{-1}$ at 21°C. The RMNS is not certified for ammonium. NO_x is derived by summing the NO₃ and NO₂ values. The assayed RMNS values per CTD deployments are listed in the appendix 8.2.

Table 10: RMNS certified concentrations \pm expanded uncertainty (U) at 21°C. Units: $\mu\text{mol L}^{-1}$

RMNS	Silicate (Si(OH) ₄)	Phosphate (PO ₄)	Nitrite (NO ₂)	NO ₃ + NO ₂ (NO _x)
Lot CP	62.5687 \pm 0.307	1.7951 \pm 0.018	0.3175 \pm 0.316	25.7136 \pm 0.379

Table 11: RMNS CP statistics for of this voyage. Units: $\mu\text{mol L}^{-1}$

RMNS CP	Silicate (Si(OH) ₄)	Phosphate (PO ₄)	Nitrite (NO ₂)	NO ₃ + NO ₂ (NO _x)
Minimum	60.800	1.780	0.313	25.700
Maximum	63.100	1.850	0.327	25.960
Mean	62.338	1.798	0.320	23.274
Median	62.500	1.790	0.319	25.830
Repeatability	0.711	0.018	0.004	0.073

7.5 Nutrient plots of RMNS

The green, pink, and red contours are at 1%, 2% and 3% from the RMNS certified mean value.

Exception: nitrite, the contours are at 0.02 $\mu\text{mol L}^{-1}$ increments from the certified value. The blue line

is the certified value's expanded uncertainty. Plots are RMNS value versus instrument run number. Please note that plots for lot CL are representing only certified values and expanded uncertainty.




Figure 7. Silicate RMNS lot CP Plot ($\mu\text{mol L}^{-1}$)




Figure 8. Phosphate RMNS lot CP Plot ($\mu\text{mol L}^{-1}$)




Figure 9. Nitrite RMNS lot CP Plot ($\mu\text{mol L}^{-1}$)




Figure 10. Nitrate + Nitrite (NOx) RMNS lot CP Plot ($\mu\text{mol L}^{-1}$)

7.6 Measurement Uncertainty

The CSIRO hydrochemistry method measurement uncertainty (MU) has been calculated for each nutrient based on the variation in the calibration curve, calibration standards, pipette and glassware calibration, and precision of the RMNS over time (Armishaw 2003).

Table 13: CSIRO Hydrochemistry nutrient analysis uncertainty values. Units: $\mu\text{mol L}^{-1}$

Calculated Measurement Uncertainty @ 1 $\mu\text{mol L}^{-1}$				
Silicate	Phosphate	Nitrite	Nitrate + Nitrite (NO _x)	Ammonia
±0.017	±0.024	±0.14	±0.019	±0.30 [‡]

The reported uncertainty is an expanded uncertainty using a coverage factor of 2 giving a 95% level of confidence.

[‡]The ammonia MU precision does not include data for the RMNS.

7.7 Method Detection Limit for Nutrients

Low nutrient seawater (LNSW) was measured 3 times in each run to determine its method detection limit (MDL). The nominal MDL was determined previously by measuring nutrients in LNSW 10 times. The MDL is set to three times the standard deviation of the LNSW results (National Association of Testing Authorities 2013). The resultant MDL was used to assess the analysis precision at low concentrations. The MDLs for each run are much lower than the nominal detection limits, indicating high analytical precision at lower concentrations. See appendix 8.3 for the measured MDL per CTD deployments.

Table 14: AA3HR auto analyser MDL statistics for this voyage. The minimum, maximum, mean, median, and reproducibility (standard deviation) are of all analytical measurements. Units: $\mu\text{mol L}^{-1}$

MDL	Silicate (Si(OH) ₄)	Phosphate (PO ₄)	Nitrate + Nitrite (NO _x)	Nitrite (NO ₂)	Ammonia (NH ₄)
Nominal MDL	0.200	0.020	0.020	0.020	0.020
Standard Dev. Min	0.30	0.09	0.01	0.01	0.16
Standard Dev. Max	0.60	0.06	0.07	0.03	0.51
Standard Dev. Mean	0.49	0.02	0.02	0.01	0.06
Standard Dev. Median	0.50	0.04	0.02	0.01	0.00
Precision of MDL (stdev)	0.09	0.05	0.02	0.01	0.22

7.8 Sampling Precision

The sampling precision for this voyage is GOOD.

Initial sampling precision is determined with the CTD test deployment (CTD 1) where multiple bottles are fired the same depth, each of which is then sampled for hydrochemistry (Table 17). Duplicate nutrient samples are also collected from the greatest depth of subsequent CTD deployments (Table 16).

For nutrients, the sampling precision is good if the difference from the mean of duplicate measurements is less than the nominal method detection limit (Table 16). The exception: NO_x (nitrate+nitrite) which uses the limit 0.06 μmol L⁻¹

Duplicate samples that exceed this limit are flagged 69 (suspect). These are tabulated in appendix 8.6.

Table 16: Difference between duplicate results. CTD 2 – CTD 5 Units: μmol L⁻¹

	Silicate (Si(OH)₄)	Phosphate (PO₄)	Nitrite (NO₂)	NO₃+ NO₂ (NO_x)	Ammonia (NH₄)
Minimum	0.200	0.000	0.001	0.050	0.000
Maximum	0.200	0.020	0.006	0.750	0.010
Mean	0.200	0.007	0.002	0.355	0.002
Variance	0.000	0.011	0.004	0.533	0.005

Table 17: CTD deployment 1. 36 bottles at 1000 dbar.

	Salinity (PSU)	Dissolved Oxygen μmol L⁻¹	Silicate (Si(OH)₄) μmol L⁻¹	Phosphate (PO₄) μmol L⁻¹	Nitrite (NO₂) μmol L⁻¹	NO₃+ NO₂ (NO_x) μmol L⁻¹
Minimum	34.439	191.551	31.200	2.040	0.001	29.550
Maximum	34.450	192.668	34.300	2.070	0.010	30.320
Mean	34.440	192.044	33.810	2.050	0.006	30.113
SD	0.002	0.282	0.620	0.005	0.003	0.157

7.9 Redfield Ratio Plot (14.0) for CTD Deployments.

The Redfield ratio for this voyage: **13.86**

The Redfield Ratio is a check for the accuracy of phosphate and nitrate+nitrite (NOx) analysis. The ratio is the required amount of P to N for marine phytoplankton growth.




Figure 13. Redfield ratio plots. Red = Data below nominal detection limit (see bad data listed in appendix 8.6).

7.10 Temperature & Humidity Change over Nutrient Analyses

The ambient conditions in the hydrochemistry laboratory and within the AA3HR instrument were measured and logged as follows:




Figure 14. Above the AA3HR instrument, temperature only. Mean 19.3°C SD 0.5°C.




Figure 15. On the deck of the AA3HR silicate and phosphate channel chemistry module, temperature only. Mean 19.6° C SD 3.3°C.

8 Appendix

8.1 Salinity: Reference material used

OSIL IAPSO Standard Seawater	
Batch	P166
Use by date	06/04/2025
K ₁₅	0.99987
PSU	34.995

8.2 Nutrients: RMNS results for each CTD Deployment.

8.2.1 Lot CP

Run analysis #	CTD Deployment #	Silicate (Si(OH) ₄) (μmol L ⁻¹)	Phosphate (PO ₄) (μmol L ⁻¹)	NOx (NO ₂ + NO ₃) (μmol L ⁻¹)	Nitrite (NO ₂) (μmol L ⁻¹)
1	1	62.550	1.790	25.863	0.315
2	2	60.933	1.783	25.743	0.325
3	3	62.567	1.793	25.862	0.320
4	4	62.500	1.793	25.770	0.326
5	5	63.067	1.837	25.750	0.316
V04_1*	NA	62.700	1.828	25.935	0.328

*In this run experimental samples were re-analysed for Phosphate.

The submitted nutrient results do **NOT** have RMNS corrections applied.

How to use the RMNS for Correction

Ratio = Certified RMNS Concentration/Measured RMNS Concentration in each run

Corrected Concentration = Ratio x Measured Nutrient Concentration

Or for smoothing data

Ratio = Average RMNS Concentration across voyage/Measured RMNS Conc. in each run

Corrected Concentration = Ratio x Measured Nutrient Concentration

8.3 Measured MDL for each CTD deployment

8.3.1 AA3HR Auto analyser

Run analysis #	CTD Deployment #	Silicate (Si(OH) ₄) (μmol L ⁻¹)	Phosphate (PO ₄) (μmol L ⁻¹)	NOx (NO ₂ + NO ₃) (μmol L ⁻¹)	Nitrite (NO ₂) (μmol L ⁻¹)	Ammonia (NH ₄) (μmol L ⁻¹)
1	1	0.600	0.043	-0.010	0.006	0.010
2	2	0.467	0.040	0.020	0.009	0.007
3	3	0.500	0.050	0.042	0.022	0.255
4	4	0.500	0.047	0.030	0.010	0.000
5	5	0.400	-0.090	0.023	0.013	-0.160
V04_1*	NA	0.125	0.000	0.040	0.010	0.010

*In this run experimental samples were re-analysed for Phosphate.

8.4 Missing or Suspect Salinity Data

Data is flagged based on CTD sampling log notes, observations during analysis, and examination of depth profile plots (Flag key: appendix 8.7)

CTD	RP	Flag	Reason for Flag
NA	NA	NA	NA

8.5 Missing or Suspect Dissolved Oxygen Data

Data is flagged based on CTD sampling log notes, observations during analysis, and examination of the depth profile (Flag key: appendix 8.7).

CTD	RP	Flag	Reason for Flag
NA	NA	NA	NA

8.6 Missing or Suspect Nutrient Data.

Data is flagged based on CTD sampling log notes, observations during analysis, and examination of the depth profile (Flag key: appendix 8.7).

CTD	RP	Analyte	Flag	Reason for Flag
1	1	Silicate	65	NA
1	3	Silicate	65	NA
1	4	Silicate	65	NA
1	5	NOx	65	NA
1	6	NOx	65	NA

8.7 Data Quality Flag Key

Flag	Description	
0	Data is GOOD	
63	Nutrients only.	Data below nominal detection limit.
65	Data is SUSPECT.	Nutrients only: Absorbance peak shape, measured by the instrument, is marginally outside set limits.
69	Data is SUSPECT.	Duplicate data is outside of set limits (software). Data point is an outlier on the depth profile plot (operator). Tagged by software or operator
79	Data is SUSPECT.	Nutrients only. Measured Method Detection Limit (MDL) for the analysis run is greater than the nominal MDL. All samples in that run tagged.
129	Data is BAD.	Nutrients Only. Absorbance peak exceeds the maximum value that can be measured by the instrument.
133	Data is BAD.	Set by operator.
134	Data is BAD.	Nutrients Only. Absorbance peak shape of calibrants, measured by the instrument, is outside of set limits (software).
141	NO Data.	Used in netcdf results file. Not used in csv results file.

8.8 GO-SHIP Specifications

8.8.1 Salinity

Accuracy of 0.001 is possible with Autosol™ salinometers and concomitant attention to methodology. Accuracy with respect to one particular batch of Standard Sea Water can be achieved at better than 0.001 PSS-78. Autosol precision is better than 0.001 PSS-78. A precision of approximately 0.0002 PSS-

78 is possible following the methods of Kawano with great care and experience. Air temperature stability of $\pm 1^\circ\text{C}$ is very important and should be recorded².

8.8.2 Dissolved Oxygen

Target accuracy is that 2 sigma should be less than 0.5% of the highest concentration found in the ocean. Precision or reproducibility (2 sigma) is 0.08% of the highest concentration found in the ocean.

8.8.3 Si(OH)_4

Approximately 1-3% accuracy¹, 0.2% precision³, full scale.

8.8.4 PO_4

Approximately 1-2% accuracy¹, 0.4% precision³, full scale.

8.8.5 NO_3

Approximately 1% accuracy¹, 0.2% precision³, full scale.

8.8.6 Notes

¹ If no absolute standards are available then accuracy should be taken to mean the reproducibility presently obtainable in the better laboratories.

² Keeping constant temperature in the room where salinities are determined greatly increases their quality. Also, room temperature during the salinity measurement should be noted for later interpretation if queries occur. Additionally, monitoring and recording the bath temperature is also recommended. The frequent use of IAPSO Standard Seawater is endorsed. To avoid the changes that occur in Standard Seawater, the use of the most recent batch is recommended. The bottles should also be used in an interleaving fashion as a consistency check within a batch and between batches.

³ Developments of reference materials for nutrients are underway that will enable improvements in the relative accuracy of measurements and clearer definition of the performance of laboratories when used appropriately and the results are reported with the appropriate meta-data.

9 References

- Armishaw, P. (2003) *"Estimating measurement uncertainty in an afternoon. A case study in the practical application of measurement uncertainty."* Accred Qual Assur, 8: pp. 218-224
- Armstrong, F.A.J., Stearns, C.A., and Strickland, J.D.H. (1967) *"The measurement of upwelling and subsequent biological processes by means of the Technicon Autoanalyzer and associated equipment,"* Deep-Sea Res., 14: pp.381-389. doi: 10.1016/0011-7471(67)90082-4
- Hood, E.M. (2010). *"Introduction to the collection of expert reports and guidelines."* The GO-SHIP Repeat Hydrography Manual: A Collection of Expert Reports and Guidelines. IOCCP Report No 14, ICPO Publication Series No. 134, Version 1, 2010.
- Hydes, D., Aoyama, M., Aminot, A., Bakker, K., Becker, S., Coverly, S., Daniel, A.G., Dickson, O., Grosso, R., Kerouel, R., van Ooijen, J., Sato, K., Tanhua, T., Woodward, E.M.S., and Zhang, J.Z. (2010). *"Determination of dissolved nutrients (N, P, Si) in seawater with high precision and inter-comparability using gas-segmented continuous flow analysers."* The GO-SHIP Repeat Hydrography Manual: A Collection of Expert Reports and Guidelines. IOCCP Report No 14, ICPO Publication Series No. 134, Version 1, 2010. (UNESCO/IOC)
- K erouel, R., and Aminot, A. (1997) *"Fluorometric determination of ammonia in sea and estuarine waters by direct segmented flow analysis"*. Mar. Chem., 57: pp. 265-275. doi: 10.1016/S0304-4203(97)00040-6
- Murphy, J. And Riley, J.P. (1962)*"A Modified Single Solution Method for the Determination of Phosphate in Natural Waters"*, Anal. Chim. Acta, 27: p.30. doi: 10.1016/S0003-2670(00)88444-5
- Rees, C., L. Pender, K. Sherrin, C. Schwanger, P. Hughes, S. Tibben, A. Marouchos, and M. Rayner. (2018) *"Methods for reproducible shipboard SFA nutrient measurement using RMNS and automated data processing."* Limnol. Oceanogr: Methods, 17(1): pp. 25-41. doi:10.1002/lom3.10294
- Wood, E.D., Armstrong, F.A.J., and Richards, F.A. (1967) *"Determination of nitrate in seawater by cadmium-copper reduction to nitrite."* Journal of the Marine Biological Association of U.K. 47: pp. 23-31.